

Characterization of neuroendocrine breast carcinomas for biomarkers of therapeutic options

Zoran Gatalica MD, DSc¹, Rebecca Feldman PhD¹, David Arguello MD¹, Anatole Ghazalpour PhD¹, Semir Vranic MD, PhD²

¹Caris Life Sciences, Phoenix, AZ; ² Department of Pathology, University of Sarajevo Clinical Center and School of Medicine, Sarajevo, Bosnia and Herzegovina



ABSTRACT

Introduction: Neuroendocrine breast carcinoma (NBC) is an extremely rare type of cancer, constituting less than 0.1% of all breast tumors, without specific treatment options. We investigated a biomarker database for frequency of molecular markers that may guide personalized treatment choices for these patients.

Materials and Methods: Molecular profiles of 40 breast carcinomas with neuroendocrine features [ER/PR+ (n=20), HER2+ (n=1) and TNBC (n=14)] were assessed (all female patients, mean age: 60.3 years, range: 39-83 years). Gene expression (Illumina DASL microarray platform), protein expression (IHC), gene amplification (ISH) and next-generation sequencing (NGS; TruSeq Illumina platform) were performed.

Results: 57% of NBCs were positive for hormone receptors (ER/PR), 40% were triple negative TNBC and 3% HER2+ subtypes. Therapeutic biomarkers (IHC) that may guide chemotherapies (and used in other primary sites neuroendocrine tumors) included: high TOP2A (85%) for etoposide or doxorubicin, low TS (57%) for 5-fluorouracil and low ERCC1 (45%) for cisplatin. Additional biomarkers for chemotherapy included: high TOPO1 (60%) for irinotecan, low RRM1 (48%) for gemcitabine and low MGMT (57%) for temozolomide. Biomarkers associated with available targeted therapies included: PTEN loss (39%), positive ALK (33%), cKit (30%), EGFR (29%), AR (26%) and PDGFRA (17%). No gene amplifications were detected in *cMET*, *EGFR*, or *TOP2A*. Targeted sequencing analysis of 47 genes detected variants in *TP53*, *PIK3CA*, *ERBB4* and *APC* genes. Gene expression data (included somatostatin receptor gene family- *SSTR1/2/3/4/5*) was available for 5 patients, for which 3/5 patients exhibited overexpression of at least one *SSTR* gene.

Conclusions: Molecular profiling by a multiplatform approach reveals potential personalized therapy options for this very rare breast cancer subtype. With recent success of somatostatin analogs for other neuroendocrine tumors, the overexpression of *SSTR* gene family in NBC is worthy of further investigation.

OBJECTIVES

Neuroendocrine breast carcinoma (NBC) is an extremely rare type of cancer, constituting less than 0.1% of all breast tumors, without specific treatment options (1). Several studies confirmed its poor prognosis (2,3).

In the present study, we investigated a biomarker database for frequency of molecular markers that may guide personalized treatment choices for these patients.


METHODS

Molecular profiles on 44 specimens were available (all female patients, mean age: 60.4 years, range: 39-83 years).

Gene expression (Illumina DASL microarray platform), protein expression (IHC), **gene amplification** (ISH) and **next-generation sequencing** (NGS; TruSeq Illumina platform) were performed.

RESULTS (UPDATED)

BC Subtypes in Cohort (n=44)



NGS genes	Mutation frequency
<i>APC</i>	7% (1/14)
<i>ERBB4</i>	7% (1/14)
<i>NOTCH1</i>	7% (1/14)
<i>PIK3CA</i>	27% (4/15)
<i>TP53</i>	50% (7/14)

Non-amplified ISH in NBC
ALK (0/1), *EGFR* (0/14), *MET* (0/8), *MYC* (0/3), *TOP2A* (0/19)

<i>SSTR</i> * mRNA expression (results)
<i>SSTR1-5</i> (1.1-3.1x)
<i>SSTR1,3,5</i> (0.8-1.3x)
<i>SSTR2,4,5</i> (0.2-12.2x)
<i>SSTR5</i> (0.4x)
<i>SSTR1,2,4</i> (0.9-2.0x)


In all, 60% (3/5) NBC exhibited gene overexpression in at least one *SSTR* (1-5) gene.

Biomarkers	Platforms	% Overexpression or amplification
ER	IHC◆	61% (27/44)
PR	IHC◆	59% (26/44)
Her-2/neu	IHC◆	2% (1/44)
<i>HER2 (ERBB2)</i>	ISH◆	0% (0/20)

◆IHC = immunohistochemistry; ISH – in-situ hybridization

* *SSTR* = Somatostatin receptor

CONCLUSIONS



* Biomarkers where no expression may predict a potential benefit.

Agents mentioned in red have favorable biomarker evidence in breast carcinoma.

REFERENCES

- Jurcic P, Kruslin B, Gatalica Z, Sanati S, Vranic S. "Breast carcinoma with neuroendocrine features: a brief review". *Endocr oncol metab* 2016;2:138-45.
- Wei B, Ding T, Xing Y, Wei W, Tian Z, Tang F, Abraham S, Nayeemuddin K, Hunt K, Wu Y. "Invasive neuroendocrine carcinoma of the breast: a distinctive subtype of aggressive mammary carcinoma". *Cancer* 2010;116:4463-73.
- Wang J, Wei B, Albarracin CT, Hu J, Abraham SC, Wu Y. "Invasive neuroendocrine carcinoma of the breast: a population-based study from the surveillance, epidemiology and end results (SEER) database". *BMC Cancer* 2014;14:147.

- Molecular profiling by a multiplatform approach reveals potential personalized therapy options for this very rare breast cancer subtype.
- Gene overexpression of the *SSTR* family in NBC is worthy of further investigation given the recent success of **somatostatin** analogs in other neuroendocrine tumors.
- The presence of ***PIK3CA*** mutations and/or loss of ***PTEN*** along with **androgen receptor** overexpression in a subset of NBCs indicates a potential use for targeted therapy (e.g. mTOR inhibitors, anti-androgen therapy) in patients with advanced and/or metastatic disease. In addition, we also identified biomarkers (***TOPO1***, ***TOPO2A***, ***RRM1***, ***ERCC1***, ***MGMT***, ***TS***) that may predict response to conventional chemotherapy options.
- Prospective clinical trials are urgently needed in this disease.